


Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GOAT ANTI-DYNAMIN 3 ANTIBODY

SKU: EB08422



SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size

Storage

Aliquot and store at -20°C. Minimize freezing and thawing. Instructions

Synonym /

Alias

dynamin 3|Dyna III|KIAA0820|MGC70433|dynamin family member

Names

Accession

NP_056384.2

Blocking

ID

Peptide

EBP08422

Immunogen

Peptide with sequence C-DKSVGNNKAEND, from the internal region of the protein sequence according to

NP_056384.2.

Peptide

Sequence

C-DKSVGNNKAEND

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography

Method using the immunizing peptide.

Shipping Instructions

Refrigerated





Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

Predicted

Human

Species Reactive

Species

Human

Human

26052 **Gene ID**

Product Grade

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

IHC Results Paraffin embedded Human Adrenal Gland. Recommended concentration: 2.5µg/ml.

ELISA

Detection

Antibody detection limit dilution 1:1000.

Limit

Blot

Western

Preliminary experiments gave an approx 26kDa band in Human Brain (Amygdala and Cerebellum), Kidney and Liver lysates after 0.3µg/ml antibody staining. Please note that currently we cannot find an explanation in the literature for the band we observe given the calculated size of 97.1kDa according to NP_056384.2. The 26kDa

band was successfully blocked by incubation with the immunizing peptide.

Application

Pep-ELISA, IHC Type

GALLERY IMAGES

