Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GOAT ANTI-EAAC1/ SLC1A1 / EAAT3 ANTIBODY

SKU: EB07181

SPECIFICATIONS

Unit Size 100 μg

Synonym / glutamate transporter|solute carrier family 1, member 1|excitatory amino acid transporter-3|excitatory amino **Alias** acid carrier 1|EAAT3|EAAC1|solute carrier family 1 (neuronal/epithelial high affinity glutamate transporter,

Names system Xag), member 1|SLC1A1

Accession

NP 004161.4 ID

Peptide with sequence C-DNEDSDTKKSYVN, from the internal region of the protein sequence according to **Immunogen**

NP 004161.4.

Peptide Sequence

C-DNEDSDTKKSYVN

Shipping

Refrigerated Instructions

Predicted

Human, Mouse, Rat, Dog **Species**

Reactive

Human **Species**

Human

6505 **Gene ID**

Mouse

20510 **Gene ID**

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png **Grade**

ELISA

Detection Antibody detection limit dilution 1:64000.

Limit

Western Western Blot: Approx 55kDa band observed in Human Placenta lysates (calculated MW of 57.1kDa according

Blot to NP_004161.4). Recommended concentration: 1-3µg/ml.

Application

Pep-ELISA, WB **Type**