


GOAT ANTI-EDD1 / HYD ANTIBODY

SKU: EB05920



SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size 100 µg

Storage Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym / Alias Names UBR5|EDD1|E3 ubiquitin protein ligase, HECT domain containing, 1|DD5|EDD|HYD|KIAA0896|progesterone induced protein|ubiquitin-protein ligase|E3 identified by differential display|hyperplastic discs protein homolog|ubiquitin protein ligase E3 component n-recogin 5|FLJ11310|MGC57263

Accession ID NP_056986.2

Blocking Peptide EBP05920

Immunogen Peptide with sequence C-LAIKTKNFGFV, from the C Terminus of the protein sequence according to NP_056986.2.

Peptide Sequence C-LAIKTKNFGFV

Purification Method Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Shipping Instructions Refrigerated

Predicted Species	Human, Mouse, Rat, Pig, Cow
Reactive Species	Human
Human Gene ID	51366
Mouse Gene ID	70790
Rat Gene ID	117060
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
IHC Results	In paraffin embedded Human Testis shows nuclear staining in all germ cells. Recommended concentration: 1-2µg/ml. Paraffin embedded Human Spleen. Recommended concentration: 1-2µg/ml.
ELISA	
Detection Limit	Antibody detection limit dilution 1:4000.
Application Type	Pep-ELISA, IHC

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

