



Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GOAT ANTI-EGR2 ANTIBODY

SKU: EB05199

SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size

Storage

Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

early growth response 2 (Krox-20 (Drosophila) homolog)|KROX-20, Drosophila, homolog (early growth Alias response-2)|Krox-20 homolog|KROX20|early growth response 2 (Krox-20 homolog, Drosophila)|EGR2 **Names**

Accession NP 000390 ID

Blocking

EBP05199 **Peptide**

Immunogen Peptide with sequence GGPLAPCSSRTRTP, from the C Terminus of the protein sequence according to NP 000390.

Peptide

GGPLAPCSSRTRTP Sequence

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using

Method the immunizing peptide.

Shipping

Refrigerated Instructions

Predicted

Human

Species

Human

1959 **Gene ID**

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring medium.png Grade

Western

Western blot: Weak band observed in 293 at approx 50kDa with 1 hour incubation of primary. Blot