

GOAT ANTI-EGR2 ANTIBODY

SKU: EB05199



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	early growth response 2 (Krox-20 (Drosophila) homolog) KROX-20, Drosophila, homolog (early growth response-2) Krox-20 homolog KROX20 early growth response 2 (Krox-20 homolog, Drosophila) EGR2
Accession ID	NP_000390
Blocking Peptide	EBP05199
Immunogen	Peptide with sequence GGPLAPCSSRTRTP, from the C Terminus of the protein sequence according to NP_000390.
Peptide Sequence	GGPLAPCSSRTRTP
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human
Human Gene ID	1959
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png
Western Blot	Western blot: Weak band observed in 293 at approx 50kDa with 1 hour incubation of primary.