

EBP06855 - IMMUNIZING PEPTIDE

SKU: EBP06855



SPECIFICATIONS

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Unit Size 100 µg

Storage Instructions Shipped at ambient temperature, store at -20°C.

Accession ID NP_644670.1

Immunogen Peptide with sequence C-SHRKSLPKAD, from the C Terminus of the protein sequence according to NP_644670.1.

Product Comments The peptide was used in the production of EBP06855 - Goat Anti-EHD4 Antibody. It is currently untested in blocking applications.

Peptide Sequence C-SHRKSLPKAD

Shipping Instructions Ambient

Product Grade https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png

Western Blot Preliminary experiments gave an approx 28-30kDa band in human breast and testes lysates after 0.05µg/ml antibody staining. Please note that currently we cannot find an explanation in the literature for the band we observe given the calculated size of 61.1kDa according to NP_644670.1. The 28-30kDa band was successfully blocked by incubation with the immunizing peptide. We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates? Have any further splice variants/modified forms been reported?

DOCUMENTS

- [Data Sheet](#)