

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GOAT ANTI-ELTD1 (DANIO RERIO AA298-310) ANTIBODY

SKU: EB11419

SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size

Storage

Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

EGF, latrophilin and seven transmembrane domain containing 1| Eltd1|zgc:63629 Alias

Names

Accession NP 998532.1 ID

Blocking

EBP11419 **Peptide**

Peptide with sequence C-EPRQAQKHPLSAN, from the internal region of the protein sequence according to Immunogen

NP 998532.1.

Peptide

C-EPRQAQKHPLSAN Sequence

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using

Method the immunizing peptide.

Shipping

Refrigerated Instructions

Predicted

Species

Zebrafish

Product

Grade

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png

ELISA

Detection Antibody detection limit dilution 1:32000.

Limit

Western

Not yet tested. At this stage we are dependent on researchers in the field for further characterization of this product. Therefore we cannot recommend an optimal concentration and the antibody is an aspiring product. We

Blot would appreciate any feedback from people in the field - have any results been reported with other

antibodies/lysates?

Application

Pep-ELISA **Type**