Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

## **GOAT ANTI-ELTD1 (DANIO RERIO AA393-407) ANTIBODY**

**SKU:** EB11420

## **SPECIFICATIONS**

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

**Unit Size** 

Storage

Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

EGF, latrophilin and seven transmembrane domain containing 1| Eltd1|zgc:63629 Alias

**Names** 

Accession NP 998532.1 ID

**Blocking** 

EBP11420 **Peptide** 

Peptide with sequence KHTQEIDPARDETKC, from the internal region of the protein sequence according to Immunogen

NP 998532.1.

**Peptide** 

KHTQEIDPARDETKC Sequence

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using

Method the immunizing peptide.

**Shipping** 

Refrigerated Instructions

**Predicted** 

Zebrafish **Species** 

**Product** 

Grade

**ELISA** 

**Detection** 

Antibody detection limit dilution 1:128000.

Limit

Western

Not yet tested. At this stage we are dependent on researchers in the field for further characterization of this product. Therefore we cannot recommend an optimal concentration and the antibody is an aspiring product. We

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring\_medium.png

Blot would appreciate any feedback from people in the field - have any results been reported with other

antibodies/lysates?

**Application** 

Pep-ELISA **Type**