

GOAT ANTI-EMP3 ANTIBODY

SKU: EB07249



SPECIFICATIONS

Unit Size	100 µg
Synonym /	
Alias	YMP HGNC:3335 epithelial membrane protein 3 EMP3
Names	
Accession	
ID	NP_001416.1
Blocking	
Peptide	EBP07249
Immunogen	Peptide with sequence DCTWNNDTKTWA, from the internal region of the protein sequence according to NP_001416.1.
Peptide	
Sequence	DCTWNNDTKTWA
Shipping	
Instructions	Refrigerated
Predicted	
Species	Human, Dog
Human	
Gene ID	2014
Product	
Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png
ELISA	
Detection	
Limit	Antibody detection limit dilution 1:128000.
Western	
Blot	Preliminary experiments in Human Kidney, Liver and Peripheral Blood Mononucleocyte lysates gave no specific signal but low background (at antibody concentration up to 1µg/ml). We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates?
Application	
Type	Pep-ELISA