


GOAT ANTI-ENDONEDD4 / NEDL2 ANTIBODY

SKU: EB06473



SPECIFICATIONS

Unit Size	100 µg
Synonym / Alias Names	NEDL2 DKFZp686M17164 EndoNEDD4 KIAA1301 protein NEDD4-related E3 ubiquitin ligase NEDL2 HECT, C2 and WW domain containing E3 ubiquitin protein ligase 2 HECW2
Usage Summary	Flow Cytometry: Flow cytometric analysis of A549 cells. Recommended concentration: 10ug/ml.
Accession ID	NP_001335697.1
Blocking Peptide	EBP06473
Immunogen	Peptide with sequence ASSAREHLLFVRR-C, from the N Terminus of the protein sequence according to NP_001335697.1.
Peptide Sequence	ASSAREHLLFVRR-C

Shipping Instructions	Refrigerated
Predicted Species	Human, Mouse, Rat, Cow
Reactive Species	Human
Human Gene ID	57520
Mouse Gene ID	227106
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
ELISA	
Detection Limit	Antibody detection limit dilution 1:32000.
Application Type	Pep-ELISA, FC

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

