
Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GOAT ANTI-ENDONEDD4 / NEDL2 ANTIBODY





SPECIFICATIONS

Unit Size 100 µg

Synonym /

Alias

Names Usage

NEDL2|DKFZp686M17164|EndoNEDD4|KIAA1301 protein|NEDD4-related E3 ubiquitin ligase NEDL2|HECT, C2 and WW domain containing E3 ubiquitin protein ligase 2|HECW2

Summary

Flow Cytometry: Flow cytometric analysis of A549 cells. Recommended concentration:

10ug/ml.

Accession

NP_001335697.1

Blocking

Peptide

EBP06473

Immunogen

Peptide with sequence ASSAREHLLFVRR-C, from the N Terminus of the protein sequence according to

NP 001335697.1.

Peptide Sequence

ASSAREHLLFVRR-C





Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

Shipping Refrigerated Instructions

Predicted

Human, Mouse, Rat, Cow

Species Reactive

Human

Species Human

Gene ID

57520

Mouse Gene ID

227106

Product

Grade

 $https://prod-vector-labs-pimcore-assets.s 3. us-east-1. amazonaws.com/assets/products/image/elite_medium.png$

ELISA

Detection

Antibody detection limit dilution 1:32000.

Limit

 $\textbf{Application} \\ \text{Pep-ELISA, FC}$ Type

GALLERY IMAGES

