

GOAT ANTI-EPOR ANTIBODY

SKU: EB06626

SPECIFICATIONS

Unit Size	100 µg
Synonym / Alias	erythropoietin receptor EPOR
Names	
Accession ID	NP_000112
Blocking Peptide	EBP06626
Immunogen	Peptide with sequence DHLGASLWPQVGC, from the N Terminus of the protein sequence according to NP_000112.
Product Comments	This antibody is expected to recognize all three reported isoforms.
Peptide Sequence	DHLGASLWPQVGC
Shipping Instructions	Refrigerated
Predicted Species	Human
Reactive Species	Human
Human Gene ID	2057
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:16000.
Western Blot	Western Blot: Approx 60kDa band observed in Human Liver lysates (calculated MW of 55.1kDa according to NP_000112). Recommended concentration: 0.3-1µg/ml.