


GOAT ANTI-ERCC6 (AA717-731) ANTIBODY

SKU: EB11429



250kDa

150kDa

100kDa

75kDa

50kDa

37kDa

25kDa

20kDa

15kDa

SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	ERCC6 Rad26 homolog RAD26 OTTHUMP00000019581 excision repair cross-complementing rodent repair deficiency, complementation group 6 DNA excision repair protein ERCC-6 CSB COFS1 COFS cockayne syndrome protein CSB Cockayne syndrome group B protein Cockayne syndrome B protein CKN2 ATP-dependent helicase ERCC6 ARMD5
Accession ID	NP_000115.1
Blocking Peptide	EBP11429
Immunogen	Peptide with sequence SNASPVQVKAYKCA, from the internal region of the protein sequence according to NP_000115.1.
Peptide Sequence	SNASPVQVKAYKCA
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human, Mouse, Rat, Dog, Pig, Cow
Reactive Species	Human
Human Gene ID	2074
Mouse Gene ID	319955
Rat Gene ID	306274
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:128000.
Western Blot	Approx 170kDa band observed in nuclear lysates of cell line HeLa (calculated MW of 168kDa according to NP_000115.1). Recommended concentration: 0.3-1µg/ml.
Application Type	Pep-ELISA, WB

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

