


GOAT ANTI-ETEA / UBXD8 ANTIBODY

SKU: EB06051



250kDa
150kDa
100kDa
75kDa
50kDa
37kDa
25kDa
20kDa
15kDa

SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	UBX domain protein 3B UBXN3B Fas associated factor family member 2 UBX domain containing 8 expressed in T-cells and eosinophils in atopic dermatitis KIAA0887 UBXD8 ETEA FAF2
Accession ID	NP_055428.1
Blocking Peptide	EBP06051
Immunogen	Peptide with sequence C-HTEVLVQDLTDE, from the C Terminus of the protein sequence according to NP_055428.1.
Peptide Sequence	C-HTEVLVQDLTDE
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human, Mouse, Rat, Dog, Cow
Reactive Species	Human, Mouse
Human Gene ID	23197
Mouse Gene ID	76577
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:16000.
Western Blot	Approx 55kDa band observed in lysates of Human and Mouse Thymus and of cell line Jurkat (calculated MW of 52.6kDa according to NP_055428.1). Recommended concentration: 0.3-1µg/ml.
Application Type	Pep-ELISA, WB

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

