Telephone: (650) 697-3600



Email: customerservice@vectorlabs.com

GOAT ANTI-F11R / JAM-A ANTIBODY

SKU: EB09489



Telephone: (650) 697-3600



250kDa 150kDa 100kDa 75kDa

50kDa

37kDa

25kDa 20kDa

15kDa

Telephone: (650) 697-3600







SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size 100 μg

Storage

Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

platelet adhesion molecule|platelet F11 receptor|junctional adhesion molecule A|junctional adhesion molecule 1|OTTHUMP00000027881|OTTHUMP00000027880|PAM-1|KAT|JCAM|JAMA|JAM1|JAM-A|JAM-1|JAM|CD321|F11

Alias Names receptor|F11R

Accession

NP_058642.1

Blocking

ID

EBP09489 **Peptide**

Immunogen

Peptide with sequence C-QPSARSEGEFKQTS, from the C Terminus of the protein sequence according to

NP 058642.1.

Peptide

C-OPSARSEGEFKOTS Sequence

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography

Method using the immunizing peptide.

Shipping

Refrigerated Instructions

Predicted

Species

Human, Mouse, Rat, Dog, Pig

Reactive

Species

Human

16456

Human

50848 Gene ID

Mouse

Gene ID

Rat Gene ID 116479

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

Grade ELISA

Detection

Western

Antibody detection limit dilution 1:2000.

Limit

Blot

Approx 38kDa band observed in Human Peripheral Blood Mononucleocyte lysates (calculated MW of 32.6kDa according to NP_058642.1). This molecular weight is routinely observed by other sources. Recommended

concentration: 0.3-2µg/ml.

Application

Pep-ELISA, WB **Type**



Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GALLERY IMAGES

250kDa 150kDa 100kDa 75kDa

50kDa

37kDa

25kDa

20kDa

15kDa