


# GOAT ANTI-F11R / JAM-A ANTIBODY

**SKU:** EB09489



250kDa  
150kDa  
100kDa  
75kDa  
  
50kDa  
37kDa  
  
25kDa  
20kDa  
  
15kDa

## SPECIFICATIONS

|                              |                                                                                                                                                                                                                                       |
|------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Formulation</b>           | Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.                                                                                                                                       |
| <b>Unit Size</b>             | 100 µg                                                                                                                                                                                                                                |
| <b>Storage Instructions</b>  | Aliquot and store at -20°C. Minimize freezing and thawing.                                                                                                                                                                            |
| <b>Synonym / Alias Names</b> | platelet adhesion molecule platelet F11 receptor junctional adhesion molecule A junctional adhesion molecule 1 OTTHUMP00000027881 OTTHUMP00000027880 PAM-1 KAT JCAM JAMA JAM1 JAM-A JAM-1 JAM CD321 F11 receptor F11R                 |
| <b>Accession ID</b>          | NP_058642.1                                                                                                                                                                                                                           |
| <b>Blocking Peptide</b>      | EBP09489                                                                                                                                                                                                                              |
| <b>Immunogen</b>             | Peptide with sequence C-QPSARSEGEFKQTS, from the C Terminus of the protein sequence according to NP_058642.1.                                                                                                                         |
| <b>Peptide Sequence</b>      | C-QPSARSEGEFKQTS                                                                                                                                                                                                                      |
| <b>Purification Method</b>   | Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.                                                                                                 |
| <b>Shipping Instructions</b> | Refrigerated                                                                                                                                                                                                                          |
| <b>Predicted Species</b>     | Human, Mouse, Rat, Dog, Pig                                                                                                                                                                                                           |
| <b>Reactive Species</b>      | Human                                                                                                                                                                                                                                 |
| <b>Human Gene ID</b>         | 50848                                                                                                                                                                                                                                 |
| <b>Mouse Gene ID</b>         | 16456                                                                                                                                                                                                                                 |
| <b>Rat Gene ID</b>           | 116479                                                                                                                                                                                                                                |
| <b>Product Grade</b>         | <a href="https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png">https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png</a>     |
| <b>ELISA Detection Limit</b> | Antibody detection limit dilution 1:2000.                                                                                                                                                                                             |
| <b>Western Blot</b>          | Approx 38kDa band observed in Human Peripheral Blood Mononucleocyte lysates (calculated MW of 32.6kDa according to NP_058642.1). This molecular weight is routinely observed by other sources. Recommended concentration: 0.3-2µg/ml. |
| <b>Application Type</b>      | Pep-ELISA, WB                                                                                                                                                                                                                         |

## DOCUMENTS

- [Data Sheet](#)

## GALLERY IMAGES

