


GOAT ANTI-FAAH ANTIBODY

SKU: EB08851



SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size 100 µg

Storage Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

Alias MGC138146|MGC102823|FAAH-1|fatty acid amide hydrolase|FAAH

Names

Accession ID NP_001432.2

Blocking Peptide EBP08851

Immunogen Peptide with sequence C-REVERLMTPEKQ, from the C Terminus of the protein sequence according to NP_001432.2.

Peptide Sequence C-REVERLMTPEKQ

Purification Method Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Shipping Instructions Refrigerated

Predicted Species	Human, Mouse, Rat, Dog
Reactive Species	Human
Human Gene ID	2166
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
IHC Results	In paraffin embedded Human Cortex shows cytoplasm staining in neuronal cell bodies. Paraffin embedded Human Testis. Recommended concentration: 3.75µg/ml.
ELISA	
Detection Limit	Antibody detection limit dilution 1:64000.
Application Type	Pep-ELISA, IHC

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

