




Telephone: (650) 697-3600



## **GOAT ANTI-FAAH ANTIBODY**

**SKU:** EB08851



## **SPECIFICATIONS**

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

**Unit Size** 

Storage

Aliquot and store at -20°C. Minimize freezing and thawing. Instructions

Synonym /

**Alias** 

MGC138146|MGC102823|FAAH-1|fatty acid amide hydrolase|FAAH

**Names** 

ID

Accession

NP\_001432.2

**Blocking** 

**Peptide** 

EBP08851

Immunogen

Peptide with sequence C-REVERLMTPEKQ, from the C Terminus of the protein sequence according to

NP\_001432.2.

**Peptide** Sequence

C-REVERLMTPEKQ

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography

Method using the immunizing peptide.

Shipping

Refrigerated Instructions





Email: <a href="mailto:customerservice@vectorlabs.com">customerservice@vectorlabs.com</a>

Telephone: (650) 697-3600

Predicted

Human, Mouse, Rat, Dog

Species Reactive

Human

Species Human

Gene ID

2166

Product Grade

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite\_medium.png

IHC Results ...

In paraffin embedded Human Cortex shows cytoplasm staining in neuronal cell bodies. Paraffin embedded

Human Testis. Recommended concentration: 3.75µg/ml.

**ELISA** 

Detection

Antibody detection limit dilution 1:64000.

Limit

Type

**Application** 

Pep-ELISA, IHC

## **GALLERY IMAGES**

