



Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GOAT ANTI-FABP4 ANTIBODY

SKU: EB12004



Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

250kDa 150kDa

100kDa 75kDa

50kDa

37kDa

25kDa

20kDa

15kDa

10kDa





Telephone: (650) 697-3600

SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size 100 µg

Storage

Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

Alias

fatty acid-binding protein, adipocyte|fatty acid-binding protein 4|adipocyte-type fatty acid-binding

protein|adipocyte lipid-binding protein|aP2|ALBP|AFABP|A-FABP|fatty acid binding protein 4, adipocyte|FABP4 Names

Accession

ID

NP_001433.1

Blocking

Peptide

EBP12004

Peptide with sequence C-SVNGDVITIKSEST, from the internal region of the protein sequence according to **Immunogen**

NP 001433.1.

Peptide

Sequence

C-SVNGDVITIKSEST

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography

Method using the immunizing peptide.

Shipping

Instructions

Refrigerated

Predicted

Species

Human, Dog

Reactive

Species

Human

Human 2167

Gene ID

Product

Grade

ELISA

Antibody detection limit dilution 1:128000.

Detection Limit

Western

Approx 14kDa band observed in Human Adipose tissue lysates (calculated MW of 14.7kDa according to

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

NP_001433.1). Recommended concentration: 0.03-0.1µg/ml.

Application

Type

Blot

Pep-ELISA, WB

GALLERY IMAGES



Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

250kDa 150kDa

100kDa 75kDa

50kDa

37kDa

25kDa

20kDa

15kDa

10kDa