

GOAT ANTI-FAZF / TZFP / ZBTB32 ANTIBODY

SKU: EB06116

250kDa

150kDa

100kDa

75kDa

50kDa

37kDa

25kDa

20kDa

15kDa

10kDa

SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	repressor of GATA OTTHUMP00000045790 ZNF538 Rog zinc finger and BTB domain containing 32 FAXF testis zinc finger protein FAZF TZFP ZBTB32
Accession ID	NP_055198.1
Blocking Peptide	EBP06116
Immunogen	Peptide with sequence RPSTSPCCPSSSTT, from the C Terminus of the protein sequence according to NP_055198.1.
Peptide Sequence	RPSTSPCCPSSSTT
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human
Reactive Species	Human
Human Gene ID	27033
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
ELISA	
Detection Limit	Antibody detection limit dilution 1:8000.
Western Blot	A faint approx 55kDa band observed in Human Bone Marrow lysates (calculated MW of 53.0kDa according to NP_055198). Recommended concentration: 2-5µg/ml.
Application Type	Pep-ELISA, WB

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

