

GOAT ANTI-FBXL3 ANTIBODY

SKU: EB07392

250kDa

150kDa

100kDa

75kDa

50kDa

37kDa

25kDa

20kDa

15kDa

10kDa

SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	FBXL3 F-box and leucine-rich repeat protein 3 FBL3 FBL3A FBXL3A F-box and leucine-rich repeat protein 3A F-box protein Fbl3a
Accession ID	NP_036290.1
Blocking Peptide	EBP07392
Immunogen	Peptide with sequence KRGGDRSDRNSSEE-C, from the N Terminus of the protein sequence according to NP_036290.1.
Peptide Sequence	KRGGDRSDRNSSEE-C
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human, Mouse, Rat, Dog, Pig
Reactive Species	Human
Human Gene ID	26224
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
IHC Results	Paraffin embedded Human Thyroid. Recommended concentration: 5µg/ml.
ELISA Detection Limit	Antibody detection limit dilution 1:32000.
Western Blot	Approx 50-55kDa band observed in Human Adipose lysates (calculated MW of 48.7kDa according to NP_036290.1). Recommended concentration: 1-3µg/ml.
Application Type	Pep-ELISA, WB, IHC

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

