



Email: <a href="mailto:customerservice@vectorlabs.com">customerservice@vectorlabs.com</a>

Telephone: (650) 697-3600

## **GOAT ANTI-FBXL3 ANTIBODY**

**SKU:** EB07392



Email: <a href="mailto:customerservice@vectorlabs.com">customerservice@vectorlabs.com</a>

Telephone: (650) 697-3600

250kDa 150kDa 100kDa

75kDa 50kDa

37kDa

25kDa 20kDa

15kDa

10kDa





Telephone: (650) 697-3600

## **SPECIFICATIONS**

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

**Unit Size** 100 μg

Storage

Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

FBXL3|F-box and leucine-rich repeat protein 3|FBL3|FBL3A|FBXL3A|F-box and leucine-rich repeat protein 3A|F-box and leuc

Alias box protein Fbl3a Names

**Accession** 

NP\_036290.1

**Blocking** 

EBP07392 **Peptide** 

Peptide with sequence KRGGRDSDRNSSEE-C, from the N Terminus of the protein sequence according to **Immunogen** 

NP 036290.1.

**Peptide** Sequence

KRGGRDSDRNSSEE-C

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography

Method using the immunizing peptide.

**Shipping** 

Refrigerated Instructions

**Predicted** 

Human, Mouse, Rat, Dog, Pig **Species** 

Reactive

Human **Species** 

Human

26224 Gene ID

**Product** 

 $https://prod-vector-labs-pimcore-assets.s 3. us-east-1. a mazon aws.com/assets/products/image/elite\_medium.png$ Grade

IHC Results Paraffin embedded Human Thyroid. Recommended concentration: 5µg/ml.

**ELISA** 

**Detection** Antibody detection limit dilution 1:32000.

Limit

Western Approx 50-55kDa band observed in Human Adipose lysates (calculated MW of 48.7kDa according to

**Blot** NP\_036290.1). Recommended concentration: 1-3µg/ml.

**Application** 

Pep-ELISA, WB, IHC **Type** 

## **GALLERY IMAGES**



Telephone: (650) 697-3600



250kDa 150kDa 100kDa

75kDa 50kDa

37kDa

25kDa 20kDa

15kDa

10kDa

