



GOAT ANTI-FBXL4 ANTIBODY

SKU: EB06044



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	F-box protein FBL4 FBL5 FBL4 F-box and leucine-rich repeat protein 4 FBXL4
Accession ID	NP_036292.2
Blocking Peptide	EBP06044
Immunogen	Peptide with sequence C-SFPKVFIKKSFTQ, from the C Terminus of the protein sequence according to NP_036292.2.
Peptide Sequence	C-SFPKVFIKKSFTQ
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human
Human Gene ID	26235
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:8000.
Western Blot	Western Blot: Preliminary experiments gave no signal but low background in Human Heart and 293 lysates at up to 1µg/ml. We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates?
Application Type	Pep-ELISA