



Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GOAT ANTI-FBXL4 ANTIBODY

SKU: EB06044

SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size

Storage

Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

Alias F-box protein FBL4|FBL5|FBL4|F-box and leucine-rich repeat protein 4|FBXL4

Names

Accession NP 036292.2 ID

Blocking

EBP06044 **Peptide**

Immunogen Peptide with sequence C-SFPKVFIKKSFTQ, from the C Terminus of the protein sequence according to NP 036292.2.

Peptide C-SFPKVFIKKSFTQ Sequence

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using

Method the immunizing peptide.

Shipping Refrigerated

Instructions

Predicted Human

Species

Human 26235

Gene ID

Product https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring medium.png Grade

Detection Antibody detection limit dilution 1:8000.

Limit

ELISA

Western Blot: Preliminary experiments gave no signal but low background in Human Heart and 293 lysates at up Western to 1µg/ml. We would appreciate any feedback from people in the field - have any results been reported with other Blot

antibodies/lysates?

Application

Pep-ELISA **Type**