

GOAT ANTI-FBXL6 ANTIBODY

SKU: EB06046



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	PP14630 FBL6A F-box protein Fbl6 FLJ22888 FBL6 F-box and leucine-rich repeat protein 6 FBXL6
Accession ID	NP_036294.2; NP_078831.4
Blocking Peptide	EBP06046
Immunogen	Peptide with sequence CLEQLLTSPSPS, from the C Terminus of the protein sequence according to NP_036294.2; NP_078831.4.
Peptide Sequence	CLEQLLTSPSPS
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human, Mouse
Human Gene ID	26233
Mouse Gene ID	30840
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:8000.
Western Blot	No signal obtained yet but low background observed in Human Kidney and Heart extracts at up to 1µg/ml. We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates?
Application Type	Pep-ELISA

DOCUMENTS

- [Data Sheet](#)