Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GOAT ANTI-FBXL6 ANTIBODY

SKU: EB06046

SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size

Storage

Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

Alias PP14630|FBL6A|F-box protein Fbl6|FLJ22888|FBL6|F-box and leucine-rich repeat protein 6|FBXL6

Names

Accession NP 036294.2; NP 078831.4 ID

Blocking

EBP06046 **Peptide**

Peptide with sequence CLEQLLTSPSPS, from the C Terminus of the protein sequence according to NP 036294.2; Immunogen

NP 078831.4.

Peptide

CLEQLLTSPSPS Sequence

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using

Method the immunizing peptide.

Shipping

Refrigerated Instructions

Predicted

Human, Mouse **Species**

Human 26233

Gene ID

Mouse

30840

Gene ID Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png Grade

ELISA

Detection Antibody detection limit dilution 1:8000.

Limit

Blot

No signal obtained yet but low background observed in Human Kidney and Heart extracts at up to 1µg/ml. We Western

would appreciate any feedback from people in the field - have any results been reported with other

antibodies/lysates?

Application

Pep-ELISA **Type**





Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600