


GOAT ANTI-FBXO11 (AA754-765) ANTIBODY

SKU: EB12386



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	viteligo-associated protein VIT-1 vitiligo-associated protein 1 ubiquitin protein ligase E3 component n-recognin 6 protein arginine N-methyltransferase 9 F-box only protein 11 VIT1 UBR6 PRMT9 FBX11 F-box protein 11 FBXO11
Accession ID	NP_079409.3; NP_001177203.1; XP_005264629.1; XP_005264630.1; XP_005264631.1; XP_005264632.1
Blocking Peptide	EBP12386
Immunogen	Peptide with sequence C-KISSYTSYPMHD, from the internal region of the protein sequence according to NP_079409.3; NP_001177203.1; XP_005264629.1; XP_005264630.1; XP_005264631.1; XP_005264632.1.
Product Comments	This antibody is expected to recognize all reported isoforms (NP_079409.3; NP_001177203.1; XP_005264629.1; XP_005264630.1; XP_005264631.1; XP_005264632.1).
Peptide Sequence	C-KISSYTSYPMHD
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human, Mouse, Rat, Pig, Cow
Reactive Species	Human
Human Gene ID	80204
Mouse Gene ID	225055
Rat Gene ID	301674
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:16000.
Western Blot	Approx 90kDa band observed in nuclear lysates of cell line HeLa and lysates of cell line A431 (calculated MW of 94.1kDa according to NP_079409.3). Recommended concentration: 0.1-0.3µg/ml. Primary incubation 1 hour at room temperature.
Application Type	Pep-ELISA, WB

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

