

GOAT ANTI-FBXO6 ANTIBODY

SKU: EB06059



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	FBXO6 F-box only protein 6 FBG2 FBX6 F-box protein FBG2 F-box protein Fbx6 F-box protein 6 RP3-330O12.3 FBG2 FBS2
Accession ID	NP_060908.1
Blocking Peptide	EBP06059
Immunogen	Peptide with sequence C-EAAQSPYRAVVQIF, from the C Terminus of the protein sequence according to NP_060908.1.
Peptide Sequence	C-EAAQSPYRAVVQIF
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human
Human Gene ID	26270
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:8000.
Western Blot	Preliminary experiments gave no signal but low background in Human Testis and Liver lysates at up to 1µg/ml. We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates?
Application Type	Pep-ELISA

DOCUMENTS

- [Data Sheet](#)