

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

## **GOAT ANTI-FBXO6 ANTIBODY**

**SKU:** EB06059

## **SPECIFICATIONS**

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

**Unit Size** 

Storage

Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

FBXO6|F-box only protein 6|FBG2|FBX6|F-box protein FBG2|F-box protein Fbx6|F-box protein **Alias** 

6|RP3-330012.3|FBG2|FBS2 **Names** 

Accession

NP 060908.1

**Blocking** 

EBP06059 **Peptide** 

Peptide with sequence C-EAAQSPYRAVVQIF, from the C Terminus of the protein sequence according to Immunogen

NP 060908.1.

**Peptide** 

C-EAAQSPYRAVVQIF Sequence

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using

Method the immunizing peptide.

**Shipping** 

Refrigerated Instructions

**Predicted** 

Human **Species** 

Human

26270 Gene ID

**Product** 

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring medium.png

Grade **ELISA** 

**Detection** Antibody detection limit dilution 1:8000.

Limit

Preliminary experiments gave no signal but low background in Human Testis and Liver lysates at up to 1µg/ml. We Western

would appreciate any feedback from people in the field - have any results been reported with other

Blot antibodies/lysates?

Application

Pep-ELISA **Type**