

GOAT ANTI-FEZ2 / ZYGIN 2 ANTIBODY

SKU: EB06331



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	zygin 2 MGC117372 fasciculation and elongation protein zeta 2 pre-T/NK cell associated protein (3Cl) fasciculation and elongation protein zeta 2 (zygin II) HUM3CL FEZ2
Accession ID	NP_005093.1
Immunogen	Peptide with sequence GERRGRDGLGAPT-C, from the N Terminus of the protein sequence according to NP_005093.1.
Peptide Sequence	GERRGRDGLGAPT-C
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human
Human Gene ID	9637
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:16000.
Western Blot	Western Blot: Preliminary experiments gave no signal but low background in Human Brain and Lung extracts at up to 1µg/ml. We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates?