

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

## **GOAT ANTI-FEZ2 / ZYGIN 2 ANTIBODY**

**SKU:** EB06331

## **SPECIFICATIONS**

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

**Unit Size** 

Storage

Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym / zygin 2|MGC117372|fasciculation and elongation protein zeta 2|pre-T/NK cell associated protein (3CI)|fasciculation **Alias** 

and elongation protein zeta 2 (zygin II)|HUM3CL|FEZ2

**Names** 

Accession NP 005093.1

Peptide with sequence GERRGRDGLGAPT-C, from the N Terminus of the protein sequence according to Immunogen NP\_005093.1.

**Peptide GERRGRDGLGAPT-C** Sequence

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using

Method the immunizing peptide.

Shipping

Refrigerated Instructions

**Predicted** 

Human **Species** 

Human

Gene ID

9637

**Product** 

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring\_medium.png Grade

**ELISA** 

**Detection** Antibody detection limit dilution 1:16000.

Limit

Western Blot: Preliminary experiments gave no signal but low background in Human Brain and Lung extracts at up Western to 1µg/ml. We would appreciate any feedback from people in the field - have any results been reported with other **Blot** 

antibodies/lysates?