

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GOAT ANTI-FGF22 ANTIBODY

SKU: EB12195

SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size

Storage

Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

Alias FGF-22|fibroblast growth factor 22|FGF22

Names

Accession NP 065688.1 ID

Blocking

EBP12195 **Peptide**

Peptide with sequence C-RVQGTRWRHGQDS, from the internal region of the protein sequence according to

Peptide

C-RVQGTRWRHGQDS Sequence

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using

Method the immunizing peptide.

Shipping

Refrigerated Instructions

Predicted

Human, Mouse, Rat **Species**

Human Gene ID

27006

Mouse **Gene ID**

67112

Rat Gene ID 170579

Product Grade

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png

ELISA

Detection Antibody detection limit dilution 1:32000.

Limit

Application Pep-ELISA

Type