

GOAT ANTI-FGF22 ANTIBODY

SKU: EB12195



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	FGF-22 fibroblast growth factor 22 FGF22
Accession ID	NP_065688.1
Blocking Peptide	EBP12195
Immunogen	Peptide with sequence C-RVQGTRWRHGQDS, from the internal region of the protein sequence according to NP_065688.1.
Peptide Sequence	C-RVQGTRWRHGQDS
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human, Mouse, Rat
Human Gene ID	27006
Mouse Gene ID	67112
Rat Gene ID	170579
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:32000.
Application Type	Pep-ELISA

DOCUMENTS

- [Data Sheet](#)