

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GOAT ANTI-FGL1 / HEPASSOCIN ANTIBODY

SKU: EB09652



Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

250kDa 150kDa

100kDa

75kDa


50kDa

37kDa

25kDa

20kDa

15kDa





Telephone: (650) 697-3600

SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size 100 μg

Storage

Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

hepatocyte-derived fibrinogen-related protein 1|hepatocellular carcinoma-related

Alias Names

sequence|hepassocin|OTTHUMP00000122468|MGC12455|LFIRE1|HP-041|HFREP1|fibrinogen-like 1|FGL1

Accession

ID

NP 004458.3

Blocking

EBP09652 **Peptide**

Peptide with sequence C-ETRVKQQQVKIKQ, from the internal region of the protein sequence according to **Immunogen**

NP 004458.3.

Product Comments

Reported variants represent identical protein: NP_963846.1, NP_963847.1, NP_671736.2, NP_004458.3

Peptide

C-ETRVKQQQVKIKQ Sequence

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography

using the immunizing peptide.

Method Shipping

Instructions

Refrigerated

Human

2267

Predicted

Species

Human, Mouse, Rat, Dog, Cow

Reactive

Species

Human

Gene ID

Gene ID

Mouse 234199

Rat Gene ID 246186

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

Grade **ELISA**

Detection

Antibody detection limit dilution 1:16000.

Limit

Western Approx 35kDa band observed in Human Liver lysates (calculated MW of 36.4kDa according to NP_004458.3).

Blot Recommended concentration: 0.1-0.3µg/ml.

Application

Type

Pep-ELISA, WB



Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GALLERY IMAGES

250kDa 150kDa

100kDa

75kDa

50kDa

37kDa

25kDa

20kDa

15kDa