

GOAT ANTI-FGL1 / HEPASSOCIN ANTIBODY

SKU: EB09652

250kDa

150kDa

100kDa

75kDa

50kDa

37kDa

25kDa

20kDa

15kDa

SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size 100 µg

Storage Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym / Alias Names hepatocyte-derived fibrinogen-related protein 1|hepatocellular carcinoma-related sequence|hepassocin|OTTHUMP00000122468|MGC12455|LFIRE1|HP-041|HFREP1|fibrinogen-like 1|FGL1

Accession ID NP_004458.3

Blocking Peptide EBP09652

Immunogen Peptide with sequence C-ETRVKQQQVKIKQ, from the internal region of the protein sequence according to NP_004458.3.

Product Comments Reported variants represent identical protein: NP_963846.1, NP_963847.1, NP_671736.2, NP_004458.3

Peptide Sequence C-ETRVKQQQVKIKQ

Purification Method Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Shipping Instructions Refrigerated

Predicted Species Human, Mouse, Rat, Dog, Cow

Reactive Species Human

Human Gene ID 2267

Mouse Gene ID 234199

Rat Gene ID 246186

Product Grade https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

ELISA Detection Limit Antibody detection limit dilution 1:16000.

Western Blot Approx 35kDa band observed in Human Liver lysates (calculated MW of 36.4kDa according to NP_004458.3). Recommended concentration: 0.1-0.3µg/ml.

Application Type Pep-ELISA, WB

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

