

GOAT ANTI-FIBRONECTIN ANTIBODY

SKU: EB11913



SPECIFICATIONS

Unit Size	100 µg
Synonym / Alias	mitigation stimulating factor(fibronectin cold insoluble globulin)OTTHUMP0000206767 OTTHUMP0000206766 OTTHUMP0000206765 OTTHUMP0000206763 OTTHUMP0000206762 OTTHUMP0000206760 OTTHUMP0000164056 OTTHUMP0000164055 OTTHUMP0000164054 OTTHUMP0000164053 H5F LTS GFND2 GFND FN2 FN NC ED-
Names	B DKF Zp686013149 DKF Zp68611370 DKF Zp68610342 DKF Zp68610164 C G fibronectin 1 FN1
Accession ID	NP_997647.1; NP_002017.1; NP_997643.1; NP_997641.1; NP_997639.1; NP_473375.2
Immunogen	Peptide with sequence KCERHTSVQTTSS, from the internal region (near N terminus) of the protein sequence according to NP_997647.1; NP_002017.1; NP_997643.1; NP_997641.1; NP_997639.1; NP_473375.2.
Product	This antibody is expected to recognize all reported isoforms (NP_997647.1; NP_002017.1; NP_997643.1; NP_997641.1; NP_997639.1; NP_473375.2).
Comments	
Peptide Sequence	KCERHTSVQTTSS
Shipping Instructions	Refrigerated
Predicted Species	Human
Human Gene ID	2335