

Telephone: (650) 697-3600

GOAT ANTI-FKBP8 / FKBP38 ANTIBODY

SKU: EB08563

SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size

Storage

Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

FK506-binding protein 8 (38kD)|FK506-binding protein 8|FK506 binding protein 8 (38kD)|FKBPr38|FKBP38|FK506 **Alias**

binding protein 8, 38kDa|FKBP8

Names Accession

NP 036313.3

Blocking

EBP08563 **Peptide**

Peptide with sequence C-KLEPSNKTIHAELSK, from the internal region of the protein sequence according to Immunogen

NP 036313.3.

Peptide Sequence

C-KLEPSNKTIHAELSK

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using

Method the immunizing peptide.

Shipping

Refrigerated Instructions

Predicted

Species

Human, Mouse, Rat, Dog, Cow

Human

23770 Gene ID

Mouse

14232

Gene ID Rat Gene ID 290652

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png

Grade

ELISA Detection

Antibody detection limit dilution 1:32000.

Limit

Application Pep-ELISA

Type