

GOAT ANTI-FKBP8 / FKBP38 ANTIBODY

SKU: EB08563



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	FK506-binding protein 8 (38kD) FK506-binding protein 8 FK506 binding protein 8 (38kD) FKBPr38 FKBP38 FK506 binding protein 8, 38kDa FKBP8
Accession ID	NP_036313.3
Blocking Peptide	EBP08563
Immunogen	Peptide with sequence C-KLEPSNKTIHAELSK, from the internal region of the protein sequence according to NP_036313.3.
Peptide Sequence	C-KLEPSNKTIHAELSK
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human, Mouse, Rat, Dog, Cow
Human Gene ID	23770
Mouse Gene ID	14232
Rat Gene ID	290652
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:32000.
Application Type	Pep-ELISA

DOCUMENTS

- [Data Sheet](#)