

GOAT ANTI-FLOTILLIN 2 / FLOT2 (INTERNAL) ANTIBODY

SKU: EB07743



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	Flotillin 2 (epidermal surface antigen 1) M17S1 ESA1 ESA reggie-1 flotillin 2 FLOT2
Accession ID	NP_004466.2
Blocking Peptide	EBP07743
Immunogen	Peptide with sequence C-QIYQDRDQFAK, from the internal region of the protein sequence according to NP_004466.2.
Peptide Sequence	C-QIYQDRDQFAK
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human, Mouse, Rat, Dog
Reactive Species	Human, Mouse, Rat
Human Gene ID	2319
Mouse Gene ID	14252
Rat Gene ID	83764
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:128000.
Western Blot	Western Blot: Approx 37kDa band observed in Human, Mouse and Rat Brain lysates (calculated MW of 41.7kDa according to NP_004466.1). Recommended concentration: 0.3-1µg/ml.

Application
Type Pep-ELISA, WB