

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GOAT ANTI-FLOTILLIN 2 / FLOT2 (INTERNAL) ANTIBODY

SKU: EB07743

SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size 100 µg

Storage

Aliquot and store at -20°C. Minimize freezing and thawing. Instructions

Synonym /

Alias Flotillin 2 (epidermal surface antigen 1)|M17S1|ESA1|ESA|reggie-1|flotillin 2|FLOT2

Names

Accession NP_004466.2 ID

Blocking

EBP07743

Peptide

Peptide with sequence C-QIYQDRDQFAK, from the internal region of the protein sequence according to **Immunogen**

NP 004466.2.

Peptide

C-QIYQDRDQFAK Sequence

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography

Method using the immunizing peptide.

Shipping

Refrigerated Instructions

Predicted

Human, Mouse, Rat, Dog

Species Reactive

Human, Mouse, Rat **Species**

Human

2319

Gene ID Mouse

14252 Gene ID

Rat Gene ID 83764

Product https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

Grade **ELISA**

Detection Antibody detection limit dilution 1:128000.

Limit

Western Western Blot: Approx 37kDa band observed in Human, Mouse and Rat Brain lysates (calculated MW of 41.7kDa

according to NP_004466.1). Recommended concentration: 0.3-1µg/ml. Blot





Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

Application Pep-ELISA, WB