



GOAT ANTI-FOXA3 / HNF3G ANTIBODY

SKU: EB06001



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	hepatocyte nuclear factor 3, gamma MGC10179 TCF3G HNF3G FKHH3 forkhead box A3 FOXA3
Accession ID	NP_004488.2
Immunogen	Peptide with sequence C-YQGLYSRLLNAS, from the C Terminus of the protein sequence according to NP_004488.2.
Peptide Sequence	C-YQGLYSRLLNAS
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human, Mouse, Rat, Zebrafish
Human Gene ID	3171
Mouse Gene ID	15377
Rat Gene ID	25100
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:2000.
Western Blot	Western Blot: No signal obtained yet but low background observed in Human Testis and Liver lysates at up to 1µg/ml. We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates?
Application Type	Pep-ELISA



www.everestbiotech.com

Email: customerservice@vectorlabs.com

Telephone: [\(650\) 697-3600](tel:(650)697-3600)