

GOAT ANTI-FOX E3 / ASMD ANTIBODY

SKU: EB06013



SPECIFICATIONS

Unit Size	100 µg
Synonym / Alias Names	forkhead-like 12 Forkhead, drosophila, homolog-like 12 forkhead (Drosophila)-like 12 FREAC8 FKHL12 ASMD forkhead box E3 FOX E3
Accession ID	NP_036318.1
Blocking Peptide	EBP06013
Immunogen	Peptide with sequence AGRSDMDPPAAFS-C, from the N Terminus of the protein sequence according to NP_036318.1.
Peptide Sequence	AGRSDMDPPAAFS-C
Shipping Instructions	Refrigerated
Predicted Species	Human
Human Gene ID	2301
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:4000.
Western Blot	No signal obtained yet but low background observed in A431 and Jurkat lysates at up to 1µg/ml. We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates?
Application Type	Pep-ELISA

DOCUMENTS

- [Data Sheet](#)

