Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GOAT ANTI-FOXE3 / ASMD ANTIBODY

SKU: EB06013

SPECIFICATIONS

Unit Size $100 \mu g$

Synonym /

forkhead-like 12|Forkhead, drosophila, homolog-like 12|forkhead (Drosophila)-like

Alias 12|FREAC8|FKHL12|ASMD|forkhead box E3|FOXE3 Names

Accession

ID

NP_036318.1

Blocking

Peptide

EBP06013

Immunogen

Peptide with sequence AGRSDMDPPAAFS-C, from the N Terminus of the protein sequence according to

NP 036318.1.

Peptide

Sequence

AGRSDMDPPAAFS-C

Shipping

Instructions

Refrigerated

Predicted

Species

Human

Human

Gene ID

2301

Product Grade

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png

ELISA

Detection

Antibody detection limit dilution 1:4000.

Limit

Western **Blot**

No signal obtained yet but low background observed in A431 and Jurkat lysates at up to 1µg/ml. We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates?

Application

Type

Pep-ELISA