

GOAT ANTI-FOXF1 ANTIBODY

SKU: EB06014



SPECIFICATIONS

Unit Size	100 µg
Synonym / Alias Names	forkhead-like 5 MGC105125 ACDMPV Forkhead, drosophila, homolog-like 5 forkhead-related activator 1 forkhead (Drosophila)-like 5 FREAC1 FKHL5 forkhead box F1 FOXF1
Accession ID	NP_001442.1
Blocking Peptide	EBP06014
Immunogen	Peptide with sequence DPASSGPSKAKKT-C, from the N Terminus of the protein sequence according to NP_001442.1.
Peptide Sequence	DPASSGPSKAKKT-C
Shipping Instructions	Refrigerated
Predicted Species	Human, Cow
Human Gene ID	2294
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:8000.
Western Blot	No signal obtained yet but low background observed in Human Kidney and Lung lysates at up to 2µg/ml. We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates?
Application Type	Pep-ELISA

DOCUMENTS

- [Data Sheet](#)

