

GOAT ANTI-FOXI2 (MOUSE) ANTIBODY

SKU: EB07010



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	Foxi2 [Mus musculus][forkhead box I2][MGI:3028075][B130055A05Rik]
Accession ID	NP_899016.1
Blocking Peptide	EBP07010
Immunogen	Peptide with sequence C-GHLYSRDGTE, from the C Terminus of the protein sequence according to NP_899016.1.
Peptide Sequence	C-GHLYSRDGTE
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Mouse, Rat
Mouse Gene ID	270004
Rat Gene ID	246073
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:8000.
Western Blot	Not yet tested. At this stage we are dependent on researchers in the field for further characterization of this product. Therefore we cannot recommend an optimal concentration and the product is investigative grade. We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates?
Application Type	Pep-ELISA

DOCUMENTS

- [Data Sheet](#)