Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

## **GOAT ANTI-FOXI2 (MOUSE) ANTIBODY**

**SKU:** EB07010

## **SPECIFICATIONS**

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

**Unit Size** 

Storage

Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

Alias Foxi2 [Mus musculus]|forkhead box I2|MGI:3028075|B130055A05Rik

**Names** 

Accession NP 899016.1 ID

**Blocking** 

EBP07010 **Peptide** 

Immunogen Peptide with sequence C-GHLIYSRDGTE, from the C Terminus of the protein sequence according to NP 899016.1.

**Peptide** 

C-GHLIYSRDGTE Sequence

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using

Method the immunizing peptide.

**Shipping** Instructions

Refrigerated

**Predicted** 

Mouse, Rat **Species** 

Mouse

270004 Gene ID

**Rat Gene ID** 246073

**Product** 

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring\_medium.png

Grade **ELISA** 

**Detection** Antibody detection limit dilution 1:8000.

Limit

Western

Not yet tested. At this stage we are dependent on researchers in the field for further characterization of this product. Therefore we cannot recommend an optimal concentration and the product is investigative grade. We

**Blot** would appreciate any feedback from people in the field - have any results been reported with other

antibodies/lysates?

Application

Pep-ELISA **Type** 





Email: <a href="mailto:customerservice@vectorlabs.com">customerservice@vectorlabs.com</a>

Telephone: (650) 697-3600