


GOAT ANTI-FOXI3 / FORKHEAD BOX I3 ANTIBODY

SKU: EB09743



SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size 100 µg

Storage Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

Alias forkhead box protein I3|FOXI3|MGI:3511278;|forkhead box I3

Names

Accession ID NP_001129121.1

Blocking Peptide EBP09743

Immunogen Peptide with sequence C-SPTISEASADTLQ, from the internal region of the protein sequence according to NP_001129121.1.

Peptide Sequence C-SPTISEASADTLQ

Purification Method Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Shipping Instructions Refrigerated

Predicted Species	Human
Reactive Species	Human
Human Gene ID	344167
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
IHC Results	In paraffin embedded Human Bronchus shows strong nuclear staining in the basal cells of the respiratory epithelium. Recommended concentration: 5-10µg/ml.
ELISA	
Detection Limit	Antibody detection limit dilution 1:4000.
Application Type	Pep-ELISA, IHC

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

