

GOAT ANTI-FOXL1 / FKH6 / FREAC7 ANTIBODY

SKU: EB06021



SPECIFICATIONS

Unit Size	100 µg
Synonym / Alias Names	forkhead-like 11 forkhead (Drosophila)-like 11 FREAC7 FKHL11 FKH6 forkhead box L1 FOXL1
Accession ID	NP_005241.1
Blocking Peptide	EBP06021
Immunogen	Peptide with sequence C-YFPLQVPDVLHFQ, from the C Terminus of the protein sequence according to NP_005241.1.
Peptide Sequence	C-YFPLQVPDVLHFQ
Shipping Instructions	Refrigerated
Predicted Species	Human
Human Gene ID	2300
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:16000.
Western Blot	Preliminary experiments gave an approx 90-100kDa band in lysates of Human Heart, Skeletal Muscle, Colon, Duodenum and Pancreas after 1µg/ml antibody staining. Please note that currently we cannot find an explanation in the literature for the band we observe given the calculated size of 36.5kDa according to NP_005241.1. The 90-100kDa band was successfully blocked by incubation with the immunizing peptide. We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates? Have any further splice variants/modified forms been reported?
Application Type	Pep-ELISA

DOCUMENTS

- [Data Sheet](#)