


GOAT ANTI-FOXO1 ANTIBODY

SKU: EB07733



SPECIFICATIONS

Unit Size	100 µg
Synonym / Alias Names	forkhead, Drosophila, homolog of, in rhabdomyosarcoma forkhead box O1A FOXO1A FKHR FKH1 forkhead box O1 FOXO1
Usage Summary	Flow Cytometry: Flow cytometric analysis of A431 cells. Recommended concentration: 10ug/ml.
Accession ID	NP_002006.2
Blocking Peptide	EBP07733
Immunogen	Peptide with sequence QDNKSSYGGMSQYNC, from the internal region of the protein sequence according to NP_002006.2.
Peptide Sequence	QDNKSSYGGMSQYNC

Shipping Instructions	Refrigerated
Predicted Species	Human
Reactive Species	Human
Human Gene ID	2308
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:128000.
Application Type	Pep-ELISA, FC

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

