

GOAT ANTI-FOXO4 / MLLT7 ANTIBODY

SKU: EB07734

250kDa

150kDa

100kDa

75kDa

50kDa

37kDa

25kDa

20kDa

15kDa

SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	myeloid/lymphoid or mixed-lineage leukemia (trithorax (Drosophila) homolog) translocated to, 7 MGC120490 AFX1 AFX myeloid/lymphoid or mixed-lineage leukemia (trithorax homolog, Drosophila) translocated to, 7 FOXO4 MLLT7
Accession ID	NP_005929.2
Blocking Peptide	EBP07734
Immunogen	Peptide with sequence C-EAASQDRMPQDLD, from the internal region of the protein sequence according to NP_005929.2.
Peptide Sequence	C-EAASQDRMPQDLD
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human
Reactive Species	Human
Human Gene ID	4303
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
ELISA	
Detection Limit	Antibody detection limit dilution 1:64000.
Western Blot	Approx 60-65kDa band observed in lysates of cell line K562 (calculated MW of 53.7kDa according to NP_005929.2). Recommended concentration: 1-3µg/ml.
Application Type	Pep-ELISA, WB

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

